

Anti-Ubiquitin antibody (STJA0040989)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB
Host / Source	Rabbit
Reactivity	Bovine/Human/Other Mammals/Non-Human Primates/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgG
Dilution Range	IHC 1:150-1:300
Formulation	Lyophilized from PBS, pH 7.2, 100mM NaCl, 15mM sodium azide.
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute in 250 μ L sterile-filtered 1X PBS, pH 7.2-7.6. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for higher stability or at 2-8°C with an appropriate antibacterial a

TARGET INFORMATION

Immunogen Sequence	Bovine
Specificity	The specificity of this antibody has been determined by indirect ELISA against ubiquitin conjugated to keyhole limpet haemocyanin. The antibody does not react with kehole limpet haemocyanin alone. Specificity has also been demonstrated by WB against

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---