

Anti-Microtubule-associated protein Tau antibody (STJA0040993)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Human/Pig
PRODUCT PROPERTIES	
Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 0.3888888888888889
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for higher stability, and at -2-8°C with an appropriate antibacterial a
TARGET INFORMATION	
Gene ID	4137
Gene Symbol	MAPT
UniProt ID	TAU_HUMAN
Immunogen Sequence	Human
Specificity	Specificity was demonstrated by Western blot. 20 ng recombinant Tau is easily detected with smaller fragments likely representing degradation products from the purified protein, which has been expressed in E. Coli. This antiserum is known to react wi
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.