

Anti-Beta-synuclein antibody (Full-Length) (STJA0040994)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:1000-1:2000 IHC 1:1000-1:2000
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C. After reconstitution divide into aliquots and store at -20°C for long-term

TARGET INFORMATION

Gene ID	6620
Gene Symbol	SNCB
UniProt ID	SYUB_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Human
Specificity	Validated by Western blotting and immunohistochemical procedures. Specific for Beta-synuclein, does not cross-react with Alpha-or Gamma-synuclein.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---