

Anti-Tyrosine Hydroxylase; Tyrosine 3-monoxygenase antibody (Full-Length) (STJA0040997)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ICC/IHC-F/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:5000-1:10000 IHC 1:100-1:500 ICC 1:100-1:500
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution of lyophilized antibody, aliquot and store at -20°C for a higher stability. Avoid freeze-thaw.

TARGET INFORMATION

Gene ID	7054
Gene Symbol	TH
UniProt ID	TY3H_HUMAN
Immunogen Region	Full-Length
Immunogen Sequence	Human
Specificity	Human. Species cross-reactivity includes rat and mouse.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---