

Anti-Neuronal nuclei antigen/RNA binding protein fox-1 homolog 3 antibody (KLH) (STJA0040999)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/ICC/IHC/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1:500-1:1000 IHC 1:5000-1:10000 ICC 1:5000-1:10000
Formulation	Lyophilized from PBS buffer pH 7.2-7.6 with 0.1% trehalose, and sodium azide
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. Store lyophilized antibody at 2-8°C After reconstitution of lyophilized antibody, aliquot and store at -20°C for

TARGET INFORMATION

Gene ID	146713
Gene Symbol	RBFOX3
UniProt ID	RFOX3_HUMAN
Immunogen Region	KLH
Immunogen Sequence	Human
Specificity	Species cross-reactivity includes human, rat, and mouse.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---