

Anti-Beta Amyloid antibody (BSA) (STJA0041000)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. At least 12 months after purchase at 2-8°C (lyophilized formulations). After reconstitution, aliquot and store at

TARGET INFORMATION

Gene ID	351
Gene Symbol	APP
UniProt ID	A4_HUMAN
Immunogen Region	BSA
Immunogen Sequence	Human
Specificity	Human beta amyloid, cross reactivity to APP has not been tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------