

## Anti-Basic Fibroblast growth factor antibody (STJA0041001)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Sheep
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:1000-1:2000 IHC 1:1000-1:2000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 $\mu$ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

### TARGET INFORMATION

Gene ID	<a href="#">2247</a>
Gene Symbol	<a href="#">FGF2</a>
UniProt ID	<a href="#">FGF2_HUMAN</a>
Immunogen Sequence	Human
Specificity	A high level of specificity for bFGF was shown by immunohistochemistry for this antiserum. This antibody is known to react with human, mouse and rat basic FGF.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---