

Anti-Neurturin antibody (STJA0041002)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC-F/WB
Host / Source	Sheep
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:2000-1:3000 IHC 1:2000-1:3000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	4902
Gene Symbol	NRTN
UniProt ID	NRTN_HUMAN
Immunogen Sequence	Human
Specificity	Dot blot shows no cross reactivity to GDNF. This antibody is known to react with human, mouse and rat. Not yet tested against other species.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---