

Anti-Pan-synuclein antibody (AKEGVVAAAEKTKQGV) (STJA0041013)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Sheep
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1 µg/mL IHC 3 µg/mL
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 µL sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	6622
Gene Symbol	SNCA
UniProt ID	SYUA_HUMAN
Immunogen Region	AKEGVVAAAEKTKQGV
Immunogen Sequence	Human
Specificity	Overlap specific immunohistochemical staining of alpha-, beta-and gamma synucleins This antiserum recognizes human and rat alpha-, beta-and gamma synucleins.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---