

Anti-Neurotrophin-3 antibody (STJA0041020)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC-F/Neut/WB
Host / Source	Sheep
Reactivity	Chicken/Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:200-1:2000 IHC 1:200-1:2000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	4908
Gene Symbol	NTE3
UniProt ID	NTE3_HUMAN
Immunogen Sequence	Human
Specificity	A cross reactivity of less than 1% to mouse NGF, recombinant human BDNF and 5% to NT4/5 has been shown by 1-site ELISA. This antiserum is known to react with rat, chicken and human NT3.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---