

Anti-Neurotrophin-3 antibody (YAEHKSHRGEY) (STJA0041024)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC-F/Neut/WB
Host / Source	Sheep
Reactivity	Chicken/Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 1-10 µg/mL IHC 1-10 µg/mL
Formulation	Lyophilized from a solution containing 3% Trehalose in PBS pH 7.4
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute in 50 µL sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	4908
Gene Symbol	NTE3
UniProt ID	NTE3_HUMAN
Immunogen Region	YAEHKSHRGEY
Immunogen Sequence	Human
Specificity	A cross reactivity of less than 0.1% to mouse NGF, recombinant human BDNF and NT4/5 has been shown by dot blot. This antiserum is known to react with rat, chicken and human NT3.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---