

Anti-Small ubiquitin-related modifier 1 antibody (aa: 6-21) (STJA0041027)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Sheep
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:4000-1:8000 IHC 1:2000-1:4000
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	7341
Gene Symbol	SUMO1
UniProt ID	SUMO1_HUMAN
Immunogen Region	aa: 6-21
Immunogen Sequence	Human
Specificity	This antiserum recognizes human SUMO-1 and not ubiquitin. This antiserum is known to cross react with rat and human SUMO-1.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---