

Anti-Gamma-synuclein antibody (114-127) (STJA0041028)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/WB
Host / Source	Sheep
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	WB 2-5 μ g/mL IHC 2-5 μ g/mL
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute in 20 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	6623
Gene Symbol	SNCG
UniProt ID	SYUG_HUMAN
Immunogen Region	114-127
Immunogen Sequence	Human
Specificity	Immunohistochemical/Western blot analysis indicate a high level of specificity for this antiserum for gamma synuclein. This antiserum is known to react with human and rat gamma synuclein.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---