

Anti-C-reactive protein antibody (STJA0041030)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB
Host / Source	Sheep
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Dilution Range	WB 1:10000-1:100000
Formulation	Lyophilized with 0.01% Thimerosal
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	1401
Gene Symbol	CRP
UniProt ID	CRP_HUMAN
Immunogen Sequence	Human
Specificity	It has been shown to be specific for C-reactive protein by WB. Human, other species have not yet been tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---