

Anti-Beta-synuclein antibody (aa: 108-125) (STJA0041034)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Sheep
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	IHC 1 µg/mL
Formulation	PBS with 0.02% sodium azide.
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 µL sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	6620
Gene Symbol	SNCB
UniProt ID	SYUB_HUMAN
Immunogen Region	aa: 108-125
Immunogen Sequence	Human
Specificity	Known to be specific for beta synuclein. This antiserum is known to react with beta synuclein of human, rat and other rodents.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---