

Anti-Alpha-synuclein antibody (aa: 108-120) (STJA0041035)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P
Host / Source	Sheep
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity-purified
Dilution Range	IHC 1 μ g/mL
Isotype	IgG
Storage Instruction	Spin vial briefly before opening. Reconstitute in 50 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	6622
Gene Symbol	SNCA
UniProt ID	SYUA_HUMAN
Immunogen Region	aa: 108-120
Immunogen Sequence	Human
Specificity	Immunohistochemical analysis of human and rat brain indicates a high level of specificity for this antiserum. Specificity was also confirmed by Western blot. This antiserum is known to react with human and rat alpha synuclein.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---