

Anti-pro-Brain-derived neurotrophic factor antibody (STJA0041036)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/Neut/WB
Host / Source	Sheep
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute in 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. After reconstitution keep aliquots at -20°C for a higher stability, and at 2-8°C with an appropriate antibacterial

TARGET INFORMATION

Gene ID	627
Gene Symbol	BDNF
UniProt ID	BDNF_HUMAN
Immunogen Sequence	Human
Specificity	Confirmed to react with purified human proBDNF, crossreact with mouse and rat proBDNF Cross reactivity with other species than human, mouse and rat has not yet been tested.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---