

Anti-IgG control antibody (STJA0041040)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Negative Control
Host / Source	Sheep
Reactivity	Non-Reactive

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	IgG
Formulation	Lyophilized from PBS, pH 7.4, without preservatives.
Isotype	IgG
Molecular Weight	14, 18, 28, 32 kDa (See application details.)
Storage Instruction	Spin vial briefly before opening. Reconstitute in 500 μ L sterile-filtered 1X PBS, pH 7.2-7.6. Centrifuge to remove any insoluble material. Store lyophilized sheep IgG at -20°C to -80°C protected from moisture. After reconstitution divide IgG into usefu

TARGET INFORMATION

Immunogen Sequence	Sheep
Specificity	Uncharacterized serum control.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---