

Anti-Serum control antibody (STJA0041041)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Negative Control
Host / Source	Sheep
Reactivity	Non-Reactive

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Whole serum
Formulation	Lyophilized, without preservatives.
Isotype	Mixed
Storage Instruction	Spin vial briefly before opening. Reconstitute with 100 μ L sterile-filtered, ultrapure water. Centrifuge to remove any insoluble material. Store lyophilized sheep whole serum at -20°C to -80°C protected from moisture. After reconstitution, divide whole

TARGET INFORMATION

Immunogen Sequence	Sheep
Specificity	Uncharacterized serum control.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---