

Anti-Annexin A2 antibody [M009] (STJA0041043)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:1000 ICC 1:100 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Molecular Weight	36
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION


Gene ID	302
Gene Symbol	ANXA2
UniProt ID	ANXA2_HUMAN
Immunogen Sequence	Human
Specificity	Clone M009 detects a 36 kDa* protein corresponding to the molecular mass of Annexin A2 on SDS-PAGE immunoblots of human A431, MDA-MB-231, MeWo, HUVEC, and PC3 cells. The antibody also works for ELISA, immunocytochemical labeling, and immunoprecipitation.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot anti-Annexin A2 antibody detection across eight human cell lysates STJA0041043



ICC staining Annexin A2 in human A431 cells using antibody with DyLight 594 detection STJA0041043