

Anti-Adiponectin antibody (C-Term) [M576] (STJA0041060)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:500
Formulation	PBS + 1 mg/ml BSA, 0.05% Na ₃ N and 50% glycerol
Isotype	IgG1
Molecular Weight	30
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	11450
Gene Symbol	Adipoq
UniProt ID	ADIPO_MOUSE
Immunogen Region	C-Term
Immunogen Sequence	Mouse
Specificity	This antibody detects a 30 kDa* protein on SDS-PAGE immunoblots of human and mouse serum, as well as detects full length human and mouse recombinant proteins.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot adiponectin detection in mouse white adipose tissue and human recombinant protein STJA0041060