

Anti-Arp3 antibody (C-Term) (STJA0041076)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IHC
Host / Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat/Xenopus

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 IHC 1:100 ICC 1:250
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	47
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION


Gene ID	10096
Gene Symbol	ACTB3
UniProt ID	ARP3_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	The antibody detects a 47 kDa* doublet corresponding to the molecular mass of Arp3 on SDS-PAGE immunoblots of human A431, Jurkat, and HeLa cells, and in rat PC12 cells.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot results anti-Arp3 antibody detection in A431, Jurkat, and HeLa cells. STJA0041076



ICC staining Arp3 in rat PC12 cells using antibody with DyLight 594 detection STJA0041076