

Anti-Argonaute 2 antibody (STJA0041077)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC
Host / Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat/Xenopus

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	97
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	27161
Gene Symbol	AGO2
UniProt ID	AGO2_HUMAN
Immunogen Sequence	Human
Specificity	The antibody detects a 97 kDa* protein corresponding to the apparent molecular mass of Ago2 on SDS-PAGE immunoblots of human A431 cells and mouse recombinant Ago2 full length protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis mouse recombinant Ago2 protein using rabbit polyclonal antibody STJA0041077