

Anti-Phospho-Agonate 2-Tyr-393 antibody (STJA0041079)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC
Host / Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat/Xenopus

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	97
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	27161
Gene Symbol	AGO2
UniProt ID	AGO2_HUMAN
Immunogen Sequence	Human
Specificity	This antibody was cross-absorbed to unphospho-Ago2 (Tyr-393) , before affinity purification using phospho-Ago2 (Tyr-393) peptide. The antibody detects a 97 kDa* protein corresponding to the apparent molecular mass of phosphorylated Ago2 on SDS-PAGE I

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis Ago2 phospho-specific antibody detection in EGF-treated A431 cells STJA0041079