

Anti-Phospho-Bad-Ser-112 antibody (STJA0041085)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500
Formulation	PBS + 0.05% NaN ₃
Isotype	IgG
Molecular Weight	23
Storage Instruction	Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20°C C or -80°C C, while keeping a working aliquot stored at 4°C C for short term. Avoid freeze/t

TARGET INFORMATION

Gene ID	12015
Gene Symbol	Bad
UniProt ID	BAD_MOUSE
Immunogen Sequence	Mouse
Specificity	The antibody detects a 23 kDa* doublet corresponding to the apparent molecular mass of phosphorylated Bad on SDS-PAGE immunoblots of mouse J774A.1 treated with calyculin A. This reactivity is not observed after lambda phosphatase treatment.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---