

Anti-CD73/NT5E antibody [M064] (STJA0041104)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:500 IP 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Molecular Weight	65
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	4907
Gene Symbol	NT5E
UniProt ID	5NTD_HUMAN
Immunogen Sequence	Human
Specificity	Clone M064 mouse monoclonal antibody detects a 65 kDa* protein on SDSPAGE Native* immunoblots of human H1299, H1915, and MDA-MB-231 carcinomas. This antibody weakly detects denatured CD73. The antibody works for native western blot, immunoprecipitati

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot CD73 detection in human cell lysates using antibody STJA0041104