

Anti-Crk II antibody (C-Term) [M332] (STJA0041117)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IP
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:1000 IP 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Molecular Weight	40
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	12928
Gene Symbol	Crk
UniProt ID	CRK_MOUSE
Immunogen Region	C-Term
Immunogen Sequence	Mouse
Specificity	This antibody detects a 40 kDa* protein corresponding to the molecular mass of Crk II on SDS-PAGE immunoblots of rat PC12 and human Jurkat cells.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



anti-Crk II antibody WB analysis detection in human Jurkat cell lysate STJA0041117