

Anti-CLASP1 antibody (C-Term) [KT66] (STJA0041128)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC
Host / Source	Rat
Reactivity	Chicken/Human/Mouse/Rat/Xenopus

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:500 ICC 1:100
Formulation	PBS + 0.05% NaN ₃
Isotype	IgG2a
Molecular Weight	161
Storage Instruction	Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival and stored long term at -20°C C or -80°C C, while keeping a working aliquot stored at 4°C C for short term. Avoid freeze/t

TARGET INFORMATION

Gene ID	76707
Gene Symbol	Clasp1
UniProt ID	CLAP1_MOUSE
Immunogen Region	C-Term
Immunogen Sequence	Mouse
Specificity	The antibody detects a 161 kDa* protein corresponding to the molecular mass of CLASP1 on SDS-PAGE immunoblots of human HeLa, mouse brain, and rat PC12 cells. The antibody also detects CLASP1 by immunocytochemistry in HeLa cells.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot anti-CLASP1 antibody detection in mouse brain and rat PC12 cells STJA0041128