

Anti-E-Cadherin antibody (C-Term) [M533] (STJA0041131)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:250
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Molecular Weight	120
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	12550
Gene Symbol	Cdh1
UniProt ID	CADH1_MOUSE
Immunogen Region	C-Term
Immunogen Sequence	Mouse
Specificity	The E-cadherin antibody detects a 120 kDa* protein corresponding to endogenous E-Cadherin expressed in human A431 cells.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot E-Cadherin detection in human A431 cells across four lanes STJA0041131