

Anti-Phospho-N-Cadherin-Tyr-820 antibody (STJA0041143)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:200
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	130
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION


Gene ID	1000
Gene Symbol	CDH2
UniProt ID	CADH2_HUMAN
Immunogen Sequence	Human
Specificity	This antibody was cross-adsorbed to phospho-N-Cadherin (Tyr-860) and unphosphorylated N-cadherin (Tyr-820) peptides before affinity purification using phospho-N-cadherin (Tyr-820) peptide. The purified antibody detects a 130 kDa* band corresponding t

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



anti-N-cadherin Tyr-820 antibody Western blot
human endothelial cells STJA0041143



ICC staining phosphorylated N-Cadherin in
pervanadate-treated mouse C2C12 cells
STJA0041143