

Anti-E-Cadherin antibody (C-Term) (STJA0041146)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	120
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	999
Gene Symbol	CDH1
UniProt ID	CADH1_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	In western blots, the antibody detects a 120 kDa band corresponding to E-cadherin in human A431 cells, and does not detect VE-cadherin or N-Cadherin.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot A431 cells probed with anti-E-Cadherin antibody at multiple dilutions
 STJA0041146