

## Anti-Phospho-Beta-Dystroglycan-Tyr-892 antibody [M117] (STJA0041171)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN <sub>3</sub> and 50% glycerol
Isotype	IgG1
Molecular Weight	43
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

### TARGET INFORMATION

Gene ID	<a href="#">1605</a>
Gene Symbol	<a href="#">DAG1</a>
UniProt ID	<a href="#">DAG1_HUMAN</a>
Immunogen Sequence	Human
Specificity	The antibody detects a 43 kDa* protein on SDS-PAGE immunoblots of human HepG2 or HeLa cells treated with pervanadate, but not in control cells.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---