

Anti-Phospho-mDia3-Ser-196 antibody (STJA0041183)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	130
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	56419
Gene Symbol	Diaph3
UniProt ID	DIAP3_MOUSE
Immunogen Sequence	Human
Specificity	The antibody detects a 130 kDa* protein corresponding to the apparent molecular mass of mDia3 on SDS-PAGE immunoblots of human Jurkat treated with calyculin A. This reactivity is not observed after lambda phosphatase treatment. The antibody may also

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot phospho-mDia3 Ser-196 detection in Jurkat cells using antibody STJA0041183