

Anti-Early Endosome Antigen 1 antibody (N-Term) [M347] (STJA0041193)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:1000 ICC 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Molecular Weight	180
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	8411
Gene Symbol	EEA1
UniProt ID	EEA1_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	This antibody detects a 180 kDa ⁺ protein corresponding to the apparent molecular mass of EEA1 on SDS-PAGE immunoblots of adult mouse brain. In immunocytochemistry, anti-EEA1 specifically stains endosomal vesicles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---