

Anti-EphA4 antibody (C-Term) (STJA0041203)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	95/120
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	2043
Gene Symbol	EPHA4
UniProt ID	EPHA4_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	The antibody detects 95 and 120 kDa* bands corresponding to EphA4 in Western blots of HUVEC cells and mouse brain tissue.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot anti-EphA4 antibodies in human endothelial cells with pervanadate treatment
 STJA0041203