

Anti-Girdin antibody (C-Term) [M012] (STJA0041232)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IHC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES



Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:250 IHC 1:50 ICC 1:50 IP 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Molecular Weight	250
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	55704
Gene Symbol	CCDC88A
UniProt ID	GRDN_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	Clone M012 mouse monoclonal antibody detects a 250 kDa* protein on SDS-PAGE immunoblots of human A431 cells and human brain tissue. The antibody also works for immunoprecipitation, immunohistochemistry, and immunocytochemistry.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot Girdin phosphorylation in ICC staining Girdin in human NCI-H1915 lung pervanadate-treated A431 cell lysates carcinoma cells using antibody STJA0041232