

Anti-Grb2 antibody (Full-Length) [M331] (STJA0041235)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IP
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:1000 ICC 1:50 IP 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Molecular Weight	24
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION


Gene ID	81504
Gene Symbol	Grb2
UniProt ID	GRB2_RAT
Immunogen Region	Full-Length
Immunogen Sequence	Rat
Specificity	This antibody detects a 24 kDa* protein corresponding to the molecular mass of Grb2 on SDS-PAGE immunoblots of rat PC12.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot anti-Grb2 antibody detection in rat PC12 cell lysate at three dilutions STJA0041235



ICC staining Grb2 in rat PC12 cells using mouse monoclonal antibody with Cy3 detection STJA0041235