

Anti-Integrin Beta 3 antibody (ECD) [M581] (STJA0041252)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	ELISA/WB/ICC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:500 ICC 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Molecular Weight	100
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.
TARGET INFORMATION	
Gene ID	3690
Gene Symbol	ITGB3
UniProt ID	ITB3_HUMAN
Immunogen Region	ECD
Immunogen Sequence	Human
Specificity	This antibody detects a 100kDa* protein corresponding to the molecular mass of Integrin Beta 3 on SDS-PAGE immunoblots of human platelets and endothelial cells (HUVEC).
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.