

Anti-I Kappa B Alpha antibody (C-Term) (STJA0041257)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	38
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	4792
Gene Symbol	NFKBIA
UniProt ID	IKBA_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	This antibody detects a 38 kDa* protein on SDS-PAGE immunoblots of human A431 and Jurkat cells, as well as mouse C2C12 cells.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---