

Anti-JMY antibody (C-Term) (STJA0041260)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC
Host / Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat/Xenopus

PRODUCT PROPERTIES



Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:200
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	130
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	57748
Gene Symbol	Jmy
UniProt ID	JMY_MOUSE
Immunogen Region	C-Term
Immunogen Sequence	Mouse
Specificity	The antibody detects a 130 kDa* protein corresponding to the molecular mass of JMY on SDS-PAGE immunoblots of human HeLa, Jurkat cells, rat PC12 cells, and mouse brain and heart.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot JMY protein expression using ICC staining JMY protein in chick fibroblasts using antibody in rathumanand mouse samples antibody green with F-actin red STJA0041260