

Anti-Memo antibody (N-Term) (STJA0041271)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC
Host / Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:200
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Molecular Weight	33
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION


Gene ID	51072
Gene Symbol	MEMO1
UniProt ID	MEMO1_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	The antibody detects a significant 33 kDa* protein corresponding to the molecular mass of Memo on SDS-PAGE immunoblots of mouse heart, mouse C2C12, and rabbit spleen fibroblasts.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot WB Memo expression in mouse and rabbit tissues using antibody STJA0041271



ICC staining Memo protein in rabbit spleen fibroblasts using antibody with DyLight 594 STJA0041271