

Anti-Na⁺/K⁺ ATPase Beta 3 antibody (ECD) [M025] (STJA0041286)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IP
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:1000 ICC 1:100 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG2a
Molecular Weight	40
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	483
Gene Symbol	ATP1B3
UniProt ID	AT1B3_HUMAN
Immunogen Region	ECD
Immunogen Sequence	Human
Specificity	Clone M025 mouse monoclonal antibody detects a 40 kDa* protein on SDS-PAGE Native* immunoblots of human A431, LNCaP, MeWo, MDA-MB-231, and MCF7 cells. This antibody does not detect denatured Na ⁺ /K ⁺ ATPase Beta 3 subunit. The antibody works for wester

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---