

Anti-Phospho-Paxillin-Tyr-31 antibody [M102] (STJA0041307)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Molecular Weight	72
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	5829
Gene Symbol	PXN
UniProt ID	PAXI_HUMAN
Immunogen Sequence	Human
Specificity	This antibody detects a 72kDa* protein corresponding to the molecular mass of phosphorylated paxillin on SDS-PAGE immunoblots of pervanadate treated A431 cells, but not in A431 control cells.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot A431 cells probed with anti-phospho-Paxillin Tyr-31 antibody STJA0041307