

Anti-PLC Gamma 1 antibody (N-Term) [M156] (STJA0041313)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IP
Host / Source	Mouse
Reactivity	Chicken/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:1000 ICC 1:100 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Molecular Weight	150
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	5335
Gene Symbol	PLCG1
UniProt ID	PLCG1_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	The antibody detects a 150 kDa* protein in human A431, Hct116, and Jurkat cells, as well as in mouse brain.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---