

Anti-PKC Delta antibody (N-Term) [M242] (STJA0041318)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/WB/ICC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:250 ICC 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG2b
Molecular Weight	78
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	5580
Gene Symbol	PRKCD
UniProt ID	KPCD_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	This antibody detects a 78 kDa* protein corresponding to the molecular mass of PKC Delta on SDS-PAGE immunoblots of adult mouse brain tissue lysate.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---