

Anti-PI3 Kinase, p85 antibody (C-Term) [M253] (STJA0041319)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC
Host / Source	Mouse
Reactivity	Chicken/Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:1000 ICC 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG2a
Molecular Weight	85
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	5295
Gene Symbol	PIK3R1
UniProt ID	P85A_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	Human
Specificity	This antibody detects an 85 kDa* band corresponding to the p85 subunit of PI3 Kinase in western blots of human A431 and Jurkat cells.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot human A431 cells probed with anti-PI3 Kinase p85 antibody at various dilutions STJA0041319



ICC staining PI3 Kinase p85 in A431 cells using antibody with DyLight 594 detection STJA0041319