

Anti-Phospho-Anti-Phosphoserine/threonine-pSer/Thr antibody [M380A/M380B] (STJA0041322)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IP
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:500 ICC 1:50 IP 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION


Immunogen Sequence	Human
Specificity	The antibody detects many serine or threonine phosphorylated proteins by western blot, immunocytochemistry, and ELISA.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western blot A431 cells phosphoserine-threonine detection with phosphatase treatment STJA0041322



ICC staining phosphoserine-phosphothreonine in A431 cells using antibody with Cy3 detection STJA0041322



Western blot NFATC1 phosphorylation in HEK-293T cells using antibody at 12000 dilution STJA0041322



Western blot WB rrdD in vitro kinase assay using serine-threonine antibody STJA0041322