

Anti-Phospho-Anti-Phosphotyrosine-p-Tyr antibody (STJA0041341)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/ICC/IP
Host / Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000 ICC 1:100 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION


Immunogen Sequence	Many
Specificity	This antibody detects phosphotyrosine-containing proteins on SDS-PAGE immunoblots of EGF treated A431 cells, as well as pervanadate treated A431, HeLa, Jurkat, and endothelial cells.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Bar graph anti-Phosphotyrosine antibody binding to phosphotyrosine peptides vs unphosphorylated control STJA0041341



ICC labeling phosphotyrosine in A431 cells using rabbit polyclonal antibody STJA0041341