

Anti-p38 Alpha MAP Kinase antibody (C-Term) (STJA0041347)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN3 and 50% glycerol
Isotype	IgG
Molecular Weight	38
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	26416
Gene Symbol	Mapk14
UniProt ID	MK14_MOUSE
Immunogen Region	C-Term
Immunogen Sequence	Mouse
Specificity	The antibody detects a 38 kDa* band corresponding to p38 Alpha on SDS-PAGE immunoblots of human Jurkat and K562, as well as mouse macrophage (J774A.1) cells.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis p38 MAP kinase using antibody in J774A1 mouse macrophage lysate STJA0041347