

Anti-Stat3 antibody (N-Term) [M263] (STJA0041380)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB/IP
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:1000 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Molecular Weight	79/86
Storage Instruction	Storage at -20°C is recommended, as aliquots may be taken without freeze/thawing due to presence of 50% glycerol. Stable for at least 1 year at -20°C.

TARGET INFORMATION

Gene ID	6774
Gene Symbol	STAT3
UniProt ID	STAT3_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	Human
Specificity	The antibody detects 86 and 79 kDa* bands corresponding to Stat3 isoforms on SDS-PAGE immunoblots of human A431 cells.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



WB analysis A431 cells Stat1Stat3and Stat5
 detection with and without EGF treatment
 STJA0041380